

Quality Certificate



Product Description:

600 ml Suspension Culture Flask Vent Cap

Quantity: 40 per case, 5 per bag
Catalog Number: 229530
Lot Number: 080411-105
Expiration Date: 20100411 (April 11, 2010)



Customer Focus:

At CELLTREAT Scientific Products, everything we do is focused on our customers. Our management and product development teams have years of industry experience. Our first hand experience working with cell and tissue culture professionals combined with our entrepreneurial spirit make us responsive to customer needs. We welcome special requests and work to provide the best and most personal service in the industry.

Satisfaction Guarantee:

CELLTREAT Scientific Products offers a total satisfaction guarantee so you can be confident in your purchasing decision. If for any reason you are not satisfied with the product performance or service provided by CELLTREAT Scientific Products we will either replace or issue a refund for the purchase price of your product.

CELLTREAT Scientific Product, LLC

P.O. Box 10
Groton, MA 01450
Phone: 877-345-4620
Fax: 978-345-4296

Quality Certificate

CELLTREAT products are manufactured in accordance with ISO 9001:2000 quality management systems in a class 100,000 clean room.

Component Materials:

Flask: Virgin Polystyrene, meets USP, Class VI requirements.

Cap: Virgin High Density Polyethylene, meets USP, Class VI requirements for plastic containers and closures.

Cap Membrane: 0.22 μ m Hydrophobic Filter

Sterilization Assurance Level: SAL 10^{-3}

This lot has been irradiated and dosimetrically released. Procedures for determining the dosage comply with AAMI/FDS 11137.

Cytotoxicity: All material resins are qualified using USP and/or ISO 10993 standards for cytotoxicity and have been shown to be non toxic.

Non-pyrogenic: Products have been validated to have less than 0.5 Eu/ml.

RNase/DNase Testing:

This product has been tested and is free of any detectable RNase/DNase contamination.

Quality Control Testing: Representative Samples have been selected and inspected according to MIL-STD-105D. This lot has passed all inspection criteria.

Country of Origin: China

This manufacturing lot has met CELLTREAT quality standards and has been accepted by quality assurance.

Shirley Happell
Quality Manager
CELLTREAT Scientific Products, LLC

CELLTREAT is a trademark of CELLTREAT Scientific Products, LLC.