

New Products from CELLTREAT!

NEW



CELL SCRAPERS

- Designed for harvesting cells from tissue culture flasks, dishes or bottles
- Available in two blade positions, "Scraper" (Part No. 229310) or "Lifter" (Part No. 229311)
- Softer, pliable blade improves cell collection
- Individually packaged
- Sterile (gamma irradiated)

Part No.	Product Description	Qty/Case
229310	Cell Scraper, 25cm Handle, 20mm "Scraper" Blade	100
229311	Cell Scraper, 25cm Handle, 20mm "Lifter" Blade	100

ERLENMEYER & FERNBACH FLASKS

- Constructed from optically clear, autoclavable polycarbonate with plain or baffled bottom
- 70mm opening for easy filling and dispensing
- Over-sized vent area in cap for improved gas exchange
- Unique 2L Fernbach now available!
- Sterile (gamma irradiated)

NEW



Part No.	Product Description	Qty/Case
229850	2L Erlenmeyer Flask, Vent Cap, Plain Bottom	6
229855	2L Erlenmeyer Flask, Vent Cap, Baffled Bottom	6
229860	2L Fernbach Flask, Vent Cap, Plain Bottom	6
229865	2L Fernbach Flask, Vent Cap, Baffled Bottom	6
229870	3L Fernbach Flask, Vent Cap, Plain Bottom	4
229875	3L Fernbach Flask, Vent Cap, Baffled Bottom	4

NEW



PLASTEUR™ PLASTIC PASTEUR PIPET

- Plasteur™ Plastic Pasteur Pipet offers the aspirating precision of a glass pipet and the safety benefits of plastic
- Economical alternative to glass
- Sterilized (gamma irradiated) so autoclaving is not necessary

Part No.	Product Description	Qty/Case
229270	Plasteur™ Plastic Pasteur Pipet, Individually Wrapped, Packed in Bags	200
229275	Plasteur™ Plastic Pasteur Pipet, Bulk packed in Bags	500
229279	Aspirating Bulb for use with Plasteur™ Pipet	12





BIO-REACTION TUBES

- Available in 15 and 50mL conical bottom and 50mL self-standing sizes
- Easy-to-read black graduations and a large white marking area
- Caps have an ultrasonically welded 0.22µm hydrophobic membrane to allow gas exchange
- Sterile (gamma irradiated)

Part No.	Product Description	Qty/Case
229471	15mL Bio-Reaction Tube - Bag	100
229472	15mL Bio-Reaction Tube - Foam Rack	300
229475	50mL Bio-Reaction Tube - Bag	100
229476	50mL Bio-Reaction Tube - Foam Rack	300
229479	50mL Self-Standing Bio-Reaction Tube - Bag	100
229473	Cap Only, Vented, Fits 15mL Bio-Reaction Tube	100
229477	Cap Only, Vented, Fits 50mL Bio-Reaction Tube	100



LARGE VOLUME CENTRIFUGE TUBES

- Available in 225mL volume
- Polyethylene cap with integral plug seal provides a secure closure
- Temperature range -80°C to 121°C
- Maximum RCF is 7,500g
- Sterile (gamma irradiated)

Part No.	Product Description	Qty/Case
229465	225mL Conical Bottom Centrifuge Tube	48



REAGENT RESERVOIRS

- Designed to assist in repetitive liquid transfer with multi-channel pipettors
- White color for solution visibility & V-shaped to reduce waste
- Sterile (gamma irradiated)

Part No.	Product Description	Qty/Case
229290	50mL Reagent Reservoir - Individual Pack	80
229291	50ml Reagent Reservoir - Bulk pack	200



For more information on the New CELLTREAT® Products listed here, or to request a sample or pricing, please visit www.celltreat.com, call 1.877.345.4620, or email us at service@celltreat.com.