

**NEW!**

# ali-Q<sup>®</sup> Edge

The **NEW** ali-Q<sup>®</sup> Edge: a **fully digital** aliquoting pipet controller delivering precision, performance and power

*Precision. Performance. Power.*

 **NEW digital display allows for exact volume selection**

LED indicators respond dynamically when aliquoting

Enter precise aliquot dispense volumes and see battery life at a glance

Easily toggle between standard and aliquoting dispense modes (0.5-5.0mL)



**Life size image!**

## ALI-Q<sup>®</sup> EDGE PIPET CONTROLLER

Part No.	Description	Aliquot Volume Range
2300-3005	ali-Q <sup>®</sup> Edge	0.5 - 5.0mL

- 1 Aspirate button
- 2 Dispense or aliquot with one button
- 3 Universal nosepiece accepts any size serological pipet
- 4 Removable kickstand
- 5 Lithium-ion battery charges with included USB-C cable
- 6 Self-standing design





# Perform equal, precise dispenses, every time!

Save both time and money from experimental errors



- Easy to handle and lightweight
- ali-Q® Edge operates as both a standard and aliquoting pipet controller — two units in one!
- Aliquot function works with any serological pipet

Compact design fits in your lab coat pocket!



Unmatched aliquoting speed and accuracy with just the push of a button



Wobble-not™ Serological Pipets are a perfect match!  
Part No. 4090-0025 is shown above

Lit Code: 110B

