

# The Reservoir You Need to Fill:

## New and Sustainable Options



### Reduce & Reuse

- Product core is made from reprocessed material, accounting for 50% of its composition
- Clearly marked with PS6 for recycling
- Domestically manufactured, reducing carbon footprint


### Cooler Accessory

- Reusable Reagent Reservoir Cooler is a custom-fit companion for VistaLab 100mL Reagent Reservoirs
- Manufactured with non-toxic gel encased in heavy duty PVC

### Small Footprint Reservoir




Small Footprint 50mL Reservoirs use 9% less plastic than standard 50mL Reservoirs



To support sustainability efforts, the inner core is comprised of reprocessed resin

Uniquely treated to provide maximum reagent recovery, up to 67% better than comparable reservoirs!

### The lowest residual volume in the industry



- 1 Trough-within-a-trough design to recover the "last few drops"
- 2 Pour-off spouts at every corner
- 3 Well-marked molded graduations on the inside wall (all sizes except 100mL)



### Mindful Packaging & Logistics

- Weekly pallets scheduled to minimize freight impact, reducing carbon emissions



### Organizations We Support



Part No.	Product Description	Qty/Case
3054-2012	10mL Reagent Reservoir, Bulk, Sterile	200
3054-2002	25mL Reagent Reservoir, Bulk, Sterile	200
3054-2006	50mL Reagent Reservoir, Bulk, Sterile	200
3054-2000	100mL Reagent Reservoir, Bulk, Sterile	200
3054-5000	100mL Reagent Reservoir Cooler	1
4054-1006-SF	<b>NEW</b> 50mL Small Footprint Reagent Reservoir, Bulk, Non-Sterile	100
3054-1006-SF	<b>NEW</b> 50mL Small Footprint Reagent Reservoir, Bulk, Sterile	100

\*Other packaging options also available



CELLTREAT Scientific Products, LLC is a Committed Sustainable Supplier



Contact your CELLTREAT rep for a sample, and/or to order

