

NEW 90mm Filter Membrane!

500mL Filter System

Product Improvement Features:


- Expanded 90mm diameter membrane increases flow and reduces clogs for faster filtration



**LARGER
AREA
SAME
PRICE!**

**SEE THE
DIFFERENCE!**

NEW Expanded
90mm Membrane
vs. 75mm Membrane



NEW

90mm Diameter
Membrane
Has 44% More
Filtration Area!

Part No.	Description	Membrane Material	Membrane Size	Filter Pore Size	Qty/ Case
229703	500mL Filter System	PES	90mm	0.45µm	12
229707	500mL Filter System	PES	90mm	0.22µm	12
229723	500mL Filter System	PES	90mm	0.10µm	12
229713	500mL Bottle Top	PES	90mm	0.45µm	24
229717	500mL Bottle Top	PES	90mm	0.22µm	24

STAY CONNECTED! Join our mailing list to stay informed of all the current CELLTREAT new product updates and promotional offers. Visit www.celltreat.com/contact to sign up.


To request a sample please email
samples@celltreat.com

CELLTREAT[®]
SCIENTIFIC PRODUCTS
www.celltreat.com



Phone: 877.345.4620 | Email: service@celltreat.com

NEW! GLASS FIBER FILTER DISKS



- Glass Fiber Filter Disks are an ideal prefilter to prevent clogging when filtering viscous or particulate-laden solutions with CELLTREAT Filter Systems and Bottle Top Filters
- Increase filtration speed by removing larger particles prior to sterilizing filtration
- Packaged sterile, individually wrapped

Part No.	Product Description	Diameter	Filter Pore Size	Qty/Case
230730	GF Filter Disk, Prefilter, Binder Free	50mm	1.0µm	50
230732	GF Filter Disk, Prefilter, Binder Free	70mm	1.0µm	50
230734	GF Filter Disk, Prefilter, Binder Free	90mm	1.0µm	50

PRODUCT IMPROVEMENT

Syringe Filters

Color coded for easy identification!



Membrane material and pore size printed on the syringe filter, opposite the luer-lock

To request a sample please email samples@celltreat.com

Phone: 877.345.4620

Email: service@celltreat.com

 **CELLTREAT**
SCIENTIFIC PRODUCTS
www.celltreat.com

