

NEW!

ali-Q[®] Edge

The **NEW** ali-Q[®] Edge: a **fully digital** aliquoting pipet controller delivering precision, performance and power

Precision. Performance. Power.



NEW digital display allows for exact volume selection

LED indicators respond dynamically when aliquoting

Enter precise aliquot dispense volumes and see battery life at a glance

Easily toggle between standard and aliquoting dispense modes (0.5-5.0mL)



Life size image!

ALI-Q[®] EDGE PIPET CONTROLLER

Part No.	Description	Aliquot Volume Range
2300-3005	ali-Q [®] Edge	0.5 - 5.0mL

- 1 Aspirate button
- 2 Dispense or aliquot with one button
- 3 Universal nosepiece accepts any size serological pipet
- 4 Removable kickstand
- 5 Lithium-ion battery charges with included USB-C cable
- 6 Self-standing design

CELLTREAT[®]
SCIENTIFIC PRODUCTS

VistaLab[™]
TECHNOLOGIES
A subsidiary of CELLTREAT[®] Scientific Products



877.345.4620



service@celltreat.com



www.celltreat.com




Perform equal, precise dispenses, every time!
Save both time and money from experimental errors



- Easy to handle and lightweight
- ali-Q® Edge operates as both a standard and aliquoting pipet controller — two units in one!
- Compatible with any serological pipet

Compact design fits in your lab coat pocket!



*Unmatched aliquoting speed and accuracy
with just the push of a button*



*Wobble-not™ Serological Pipets are a perfect match!
Part No. 4090-0025 is shown above*

Lit Code: 110A

