



OVATION LfTS PIPETTE TIP FIT GUIDE

Find the perfect fit tip for your NEW Ovation LfTS

Ovation - Electronic Single Channel						
Part Number	Volume Range	Filter LfTS Tip	Non-filter LfTS Tip	Filter LfTS Tip	Non-filter LfTS Tip	
1065-0020LF	0.5 - 20µL	229071	229074	4080-1332	4080-1032	
1065-0125LF-200	2 - 125µL	229072	229075	4080-2332	4080-2032	
1065-0125LF-300	2 - 125µL	229081	229083	4080-2333	4080-2033	
1065-0250LF-200	5 - 250µL	NA	229075	NA	4080-2032	
1065-0250LF-300	5 - 250µL	229081	229083	4080-2333	4080-2033	
1065-1250LF	25 - 1250µL	229082	229084	4080-3333	4080-3033	
Ovation - Electronic Multichannel						
Part Number	Volume Range	Filter LfTS Tip	Non-filter LfTS Tip	Filter LfTS Tip	Non-filter LfTS Tip	
1160-0020LF	0.5 - 20µL	229071	229074	4080-1332	4080-1032	
1160-0125LF-200	2 - 125µL	229072	229075	4080-2332	4080-2032	
1160-0125LF-300	2 - 125µL	229081	229083	4080-2333	4080-2033	
1160-0250LF-200	5 - 250µL	NA	229075	NA	4080-2032	
1160-0250LF-300	5 - 250µL	229081	229083	4080-2333	4080-2033	
1160-1250LF	25 - 1250µL	229082	229084	4080-3333	4080-3033	
1060-0020LF	0.5 - 20µL	229071	229074	4080-1332	4080-1032	
1060-0125LF-200	2 - 125µL	229072	229075	4080-2332	4080-2032	
1060-0125LF-300	2 - 125µL	229081	229083	4080-2333	4080-2033	
1060-0250LF-200	5 - 250µL	NA	229075	NA	4080-2032	
1060-0250LF-300	5 - 250µL	229081	229083	4080-2333	4080-2033	
1060-0850LF	25 - 850µL	229073	229076	4080-3332	4080-3032	

Ovation - Mechanical					
Part Number	Volume Range	Filter LfTS Tip	Non-filter LfTS Tip	Filter LfTS Tip	Non-filter LfTS Tip
1070-0002LF	0.2 - 2µL	229071	229074	4080-1332	4080-1032
1070-0010LF	1 - 10µL	229071	229074	4080-1332	4080-1032
1070-0020LF	2 - 20µL	229071	229074	4080-1332	4080-1032
1070-0100LF-200	10 - 100µL	229072	229075	4080-2332	4080-2032
1070-0100LF-300	10 - 100µL	229081	229083	4080-2333	4080-2033
1070-0200LF-200	20 - 200µL	229072	229075	4080-2332	4080-2032
1070-0200LF-300	20 - 200µL	229081	229083	4080-2333	4080-2033
1070-1000LF	100 - 1000µL	229073	229076	4080-3332	4080-3032

Ovation - Quick Set					
Part Number	Volume Range	Filter LfTS Tip	Non-filter LfTS Tip	Filter LfTS Tip	Non-filter LfTS Tip
1057-0002LF	0.2 - 2µL	229071	229074	4080-1332	4080-1032
1057-0010LF	1 - 10µL	229071	229074	4080-1332	4080-1032
1057-0020LF	2 - 20µL	229071	229074	4080-1332	4080-1032
1057-0100LF-200	10 - 100µL	229072	229075	4080-2332	4080-2032
1057-0100LF-300	10 - 100µL	229081	229083	4080-2333	4080-2033
1057-0200LF-200	20 - 200µL	229072	229075	4080-2332	4080-2032
1057-0200LF-300	20 - 200µL	229081	229083	4080-2333	4080-2033
1057-1000LF	100 - 1000µL	229073	229076	4080-3332	4080-3032

To ensure there is an Ovation that fits every lab, we offer certain models with the option of the 200µL or 300µL LfTS fit tip.

Didn't order the correct fit, or your protocols or tip size ordering has changed?

No problem! Send in your Ovation and we can swap nozzles for you. Please email TechService@celltreat.com.