



D-Tipper Pinks

- Special Edition Pinks supporting Breast Cancer Research
- Fixed volume pipettes
- Single-stroke “to deliver” operation
- Low force plunger
- Built-in tip ejector
- Glove guard
- In-lab calibration adjustment

\$5 will be donated to Memorial Sloan Kettering Cancer Center for each Pinks pipette sold.



| VistaLab Catalog # | Color | Volume | Accuracy | Precision |
|--------------------|-------|--------|----------|-----------|
| 1051PK | Pink | 10µL | ±5.00% | 2.60% |
| 1052PK | Pink | 20µL | ±3.00% | 1.35% |
| 1054PK | Pink | 50µL | ±2.00% | 0.80% |
| 1055PK | Pink | 100µL | ±1.00% | 0.54% |
| 1056PK | Pink | 200µL | ±1.00% | 0.30% |
| 1057PK | Pink | 250µL | ±1.00% | 0.72% |
| 1059PK | Pink | 500µL | ±1.00% | 0.40% |
| 1060PK | Pink | 1000µL | ±1.00% | 0.60% |